

NEW TFT Display

Winstar, introduced the new 15" IPS TFT-LCD module with a resolution of 1024x768 pixels. This TFT module adopted IPS technology which has the advantage of a wider viewing angle of Left:88 / Right:88 / Up:88 / Down:88 degrees. it supports LVDS interface, with brightness of 300 nits (typical value), contrast ratio of 2000:1, anti-glare surface glass panel.

Operational Temperature Range: -20°C to +70°C

Storage Temperature Range: -30°C to +70°C.

WF150ATYAMLNNO

Item	Contents	Unit
Size	15	inch
Dot Matrix	1024 × RGB × 768 (TFT)	dots
Module Dimension	326.5 × 253.5 × 9.1	mm
Active Area	304.1 × 228.1	mm
Dot Pitch	0.297 × 0.297	mm
LCD Type	TFT, Normally Black, Transmissive	----
Viewing Angle	88/88/88/88	----
Backlight Type	LED, Normally White	----
Interface	LVDS	----
Touch Panel	Without Touch Panel	----
Surface	Anti-Glare	----

The image shows a medical clinic status query interface on a TFT display. The title bar is green and reads "QUERY ON CLINIC STATUS". Below the title bar, there are two main sections: "Query by Clinic Department" and "Query by patient information".

Query by Clinic Department: This section includes a text input field for "Department No.", a "Please select a consultation time" section with three buttons: "Morning", "Afternoon", and "Night", and "OK" and "Cancel" buttons.

Query by patient information: This section includes a text input field for "Chart No / Passport Number", a "Date of Birth (EX: YY/MM/DD)" section with three input fields for YY, MM, and DD, and "OK" and "Cancel" buttons.

Below these sections is a numeric keypad (0-9) and an alphanumeric keypad (A-Z). The numeric keypad has a "Delete" button, and the alphanumeric keypad has a "Backspace" button.

At the bottom of the screen, there is a row of six circular icons with labels: "Registered service", "Reschedule or Cancel your appointment", "Query on Clinic Status", "Departments", "Transportation", and "Health Examination".